

Brewer Spectrophotometer laboratory wavelength characterization with a tunable-laser

- implications to the Brewer ozone retrieval-

Alberto Redondas¹, Saulius Nevas², Alberto Berjón¹, Meelis-Mait Sildoja², Virgilio Carreño¹, Daniel Santana¹, Sergio F. León-Luis¹, Bentorey Hernández¹ & Javier López-Solano¹

¹ *Izaña Atmospheric Research Center, AEMET, Tenerife, Canary Islands, Spain*

² *Physikalisch-Technische Bundesanstalt, 38116, Braunschweig, Germany*



#1 - Dispersion measurements using spectral lamps (zinc, cadmium and mercury).

#2- Dispersion measurements using tunable laser (Brewer scanning).


While Laser is emitting in a fixed wavelength, the Brewer will scan around this wavelength ($\pm 2\text{nm}$) moving the grating and using the 6 slits. This process will be repeated 16 times, using wavelengths ranging from 290 nm to 365 nm with an increment of 5 nm.



#3 - Dispersion measurements using tunable laser (Laser scanning).



While Brewer is measuring in ozone mode and in aerosol mode, the laser will scan around ($\pm 2\text{nm}$) the six Brewer slits with a step of 0.05nm for different grating positions.



- Brewer in ozone mode. Laser scan $\pm 2\text{nm}$ from the 6 brewer slits in ozone mode, with a step of 0.05 nm.
- Brewer in **aerosol mode**. Laser scan $\pm 2\text{nm}$ from the 6 brewer slits in aerosol mode (4 different grating positions), with a step of 0.05 nm.


- Setup for bandpass and wavelength characterisation of Brewer at PTB











$$O_3 = \frac{F - ETC}{\alpha m}$$

Wavelength Calibration

1. Optimized wavelengths (Sun Scan)
2. Dispersion test, (discharge lamps)
3. Ozone Cross Sections (Bass & Paur, $T = -45$).

$$\tilde{\alpha}(X, \mu) = \sum w_i \frac{\int \alpha(\lambda) * S(\lambda, \lambda') * F(\lambda, \lambda', X, \mu) d\lambda'}{\int S(\lambda, \lambda') * F(\lambda, \lambda', X, \mu) d\lambda'}$$

Alpha = ozone cross section

Wi= ozone weighting coefficients

S = Slit Function, (ILS)

F = Sun Spectra , (X, and nu)

- Use “ideal” slits , (no wings, no stray light, '0' outside the triangle)
- The FWHM of the triangle is dependent of the slit.
- This o3 cross section is Bass & Paur absorption coefficient.
- Flat Spectra : Solar Spectra is not considered
- The O3abs (also SO2, Rayleigh,) are calculated for an interval of steps around the ozone position, for possible cal-step changes.







Brewer Dispersion procedure summary:

The brewer scan every 10 steps (~ 0.6) Å the emission lines of Hg and Cd, using a particular slit

1. Central step and fwhm (step) are calculated for the six slits, assuming a isosceles triangle.
2. The dispersion relation (wavelength vs step) is determined
 - a. Standard quadratic polynomial
 - b. Cubic (Gröebner & Ker 1998)
3. The central and fwhm steps are converted to wavelength.
4. Parametrized Slit are used for cross section calculation.

Gröbner, J., Wardle, D. I., McElroy, C. T. y Kerr, J. B.: Investigation of the wavelength accuracy of Brewer spectrophotometers, APPLIED OPTICS, 37(36), 1998.

Brewer Dispersion (I)



Lamp	Line (nm)	Slits
Mercury (Hg)	289.36	0-1
Indium (In)	293.263	0-3
Hg	296.728	0-3
Zinc (Zn)	301.836	0-5
Zn	303.578	0-5
Cd (multiplet)	308.082	0-5
Cd	313.3167	0-5
Cd	326.1055	0-5
Zn	328.233	0-5
Hg	334.148	0-5
Cd	340.3652	0-5
Cd	349.995	4-5
Cd (multiplet)	361.163	5

Saving ozonecoefs to ./opos02616_PTB.185

1020 WL(A)	3031.84	3062.95	3100.43	3134.90	3167.98	3199.90	
Res(A)	11.39	11.12	10.93	11.24	10.98	10.83	
O3abs(1/cm)	2.6018	1.7813	1.0051	0.6765	0.3750	0.2943	O3: 0.3423
Daumt O3abs(1/cm)	2.6086	1.7761	0.9968	0.6706	0.3681	0.2949	O3: 0.3530
So2abs(1/cm)	3.4603	5.6249	2.4003	1.9180	1.0545	0.6144	
1e4*Rayabs(1/cm)	5051.1	4832.6	4585.2	4371.7	4178.7	4002.6	R: 0.0011
I0(mW m ⁻² nm ⁻¹)	0.07	0.12	0.19	0.32	0.40	0.48	

Ozone offset due to Rayleigh: -3.1 DU


Ratio Ozone for So2(A3)= 1.1483, So2/O3(A2)= 2.7537



Non linearity of PMT

- The PMT requires a different configuration to measure pulsed sources (thanks Volodya !!).
- We can work on “linear regions” using the direct port for maximizing the laser output.


Two critical areas:







- Low counts ~ 100 c/s
- Hysteresis $\sim 10^4$ c/s






Laser scan : brewer measure in ozone mode (six slits + dark)






Dark correction




- No linearity and dark correction is not applied
- Slit # 1 measurements correspond to the DARK (blocked) position
- In some experiments the **dark count** is higher than the measurements the brewer correction assign two counts for negative values
- The hysteresis of the no linearity is clear for the asymmetry of the peak.




- Dark values were highest immediately after exposing the PMT to the laser light. The dark signal was then gradually decreasing with time.
- The dark values of the PMT may show a fast decay after the excitation, which may cause the values for slit #1 (measured immediately after slit #0) be higher than for some other slits, measured later.


The best results are obtained if the dark count is also corrected by nonlinearity and then subtracted



But,.. only works with slit#0




Laser Scanning : Ozone and AOD mode









Brewer Slit Parametrization

The brewer algorithms assume a trapezoidal slit, cut a 0.87, center and fwhm at the calculated for every slit.

Using the measured slit, the effective ozone cross section calculated is **~0.9% higher** respect the parametrized brewer slit used on the standard procedure.



	Parametrized (B&P)	Measured (B&P)
Brw	0.3381	0.3407
B&P	0.333	0.336
DMB	0.3483	0.3514
SDK	0.3392	0.3422



Mic	Slit #	wv	fwhm	Measured	Parametrize
1020	0	3032.0271	5.7258	6.0007	5.9989
1020	2	3063.0381	5.5522	4.0936	4.0935
1020	3	3100.5519	5.4846	2.3063	2.3065
1020	4	3135.0695	5.7049	1.5652	1.565
1020	5	3168.0487	5.5874	0.8598	0.8602
1020	6	3200.0841	5.5148	0.68	0.6784

Ozone Absorption Coefficient


Matlab date	mic	slit	brw_wv	brw_fwhm	brw_sdk_me	brw_sdk_par
736355.7552	1020	0	3032.02711	5.72579	6.00073	5.9989
736355.7734	1020	2	3063.03809	5.55223	4.09356	4.0935
736355.6786	1020	3	3100.55194	5.48458	2.30626	2.30654
736355.7006	1020	4	3135.06955	5.7049	1.5652	1.56504
736355.7187	1020	5	3168.04868	5.58735	0.85981	0.86016
736355.7369	1020	6	3200.08413	5.51482	0.68	0.67842
		Ozone Ratio		0.34226	0.34091	

ozone wv	brw_meas	brw_param
	0.34226	0.34091
ratio	1.00398	1

~0.4 A in FWHM produces 1.2% in the retrieved ozone


Brewer OPO Scan

-
- Laser emits 295 to 363 every 5 nm
- Brewer scan every 5 steps resolution (~ 0.3 A) (doubled)
- The readings are corrected by nonlinearity




laser_ptb/BrewerScan_0160129.txt

362.9938 nm slit=5 line=15





laser_ptb/BrewerScan_20160129.txt

362.9938 nm slit=5 line=15




Residuals of the fit




Wavelength calibration experiments :

1. **Laser** wavelength scanned @ fixed Brewer position : **Laser_brw_o3**
2. Brewer grating pos. changed @ fixed **laser** wavelength, dispersion approx. by a **quadratic** function : **Brewer_laser_q**
3. Brewer grating pos. changed @ fixed **laser** wavelength, dispersion approx. by a **cubic** function : **Brewer_laser_c**
4. Brewer grating pos. changed @ fixed **lamp** emission wavelength, dispersion approx. by a **quadratic** function: **Brewer_lamp_q**
5. Brewer grating pos. changed @ fixed **lamp** emission wavelength, dispersion approx. by a **cubic** function : **Brewer_lamp_c**




FWHM vs wavelength


Central wavelength diff




Central wavelength diff




FWHM wavelength diff




FWHM wavelength diff



Serdyuchenco effective ozone cross section ratio



Serdyuchenco effective ozone cross section ratio



Wavelength calibration methods:

1. **Laser** wavelength scanned @ fixed Brewer position
2. Brewer grating pos. changed @ fixed **laser** wavelength, dispersion approx. by a **quadratic** function
3. Brewer grating pos. changed @ fixed **laser** wavelength, dispersion approx. by a **cubic** function
4. Brewer grating pos. changed @ fixed **lamp** emission wavelength, dispersion approx. by a **quadratic** function
5. Brewer grating pos. changed @ fixed **lamp** emission wavelength, dispersion approx. by a **cubic** function

Method	1.	2.	3.	4.	5.
SGW	0.3409	0.3442	0.342	0.3446	0.3412
ratio	1	1.0096	1.0033	1.0108	1.001

Título de la presentación

Many thanks :

Volodya Savastiouk

Pavel Babal

Keith Wilson

Still data for analyze

- Laser replication of Cd lines.
- Repeated laser measurements.
- Comparison for HeCd laser.
-

This work has been supported by the European Metrology Research Programme (EMRP) within the joint research project ENV59 “Traceability for atmospheric total column ozone” (ATMOZ). The EMRP is jointly funded by the EMRP participating countries within EURAMET and the European Union.



Presentación Alberto

The RBCC-E Team



Alberto Redondas (AEMET)
Alberto Berjon (ULL, ATMOZ)
Javier López Solano (ULL, IDEAS)
Bentorey Hernandez (ULL, PANDONIA)
Virgilio Carreño (AEMET)
Manuel Rodriguez Valido (ULL)
Daniel Santana (ULL, PANDONIA)
Sergio Fabián León Luis (AEMET)

Filiación Completa

Alberto Redondas Marrero, RBCC-E site manager Izaña Atmospheric Research Centre, State Meteorological Agency of Spain.

Virgilio Carreño Corbella , RBCC-E researcher, Izaña Atmospheric Research Centre, State Meteorological Agency of Spain.

Alberto Berjón Arroyo, RBCC-E researcher, University of La Laguna, Spain.