

Temperature characterisation of Brewer determined in the laboratory

Alberto Berjón, Alberto Redondas, Meelis-Mait Sildoja, Saulius Nevas, Virgilio Carreño, Daniel Santana, Sergio F. León-Luis, Bentorey Hernández, Javier López-Solano



This work has been supported by the European Metrology Research Programme (EMRP) within the joint research project ENV59 “Traceability for atmospheric total column ozone” (ATMOZ). The EMRP is jointly funded by the EMRP participating countries within EURAMET and the European Union.

Objectives:

- Determine the effect of temperature on the Brewer measurements made through different ports: Direct and Global ports, and internal lamp.
- Determine the effect of PTFE diffuser on the instrument temperature dependence.
- Compare coefficients obtained in the laboratory and using field measurements.

Brewer TOC calculation

Beer-Lambert law:
$$TOC = \frac{R_6 - ETC - B}{A\mu}$$

Intensity corrections:

- Dark counts
- Deadtime
- Temperature
- Filter transmittance

$$R_6 = \sum_{i=1}^n w_i F(\lambda_i)$$

$$F(\lambda_i) = 10^4 \log(I(\lambda_i))$$

$$ETC = \sum_{i=1}^n w_i F_0(\lambda_i)$$

$$B = \nu \frac{p}{p_0} \sum_{i=1}^n w_i \beta(\lambda_i)$$

$$A = \sum_{i=1}^n w_i \alpha(\lambda_i)$$

Temperature correction

$$I = I_c - \tau_0(T - T_0) \quad I_c = \frac{I}{1 - \tau(T - T_0)} \quad \tau = \tau_0/I_c$$

$$\ln(I_c) = \ln(I) + \tau(T - T_0)$$

Temperature correction

$$\ln(I_c) = \ln(I) + \tau (T - T_0)$$

$$10^4 \log(I_c) = 10^4 \log(I) + \tau_b T \quad F_c = F + \tau_b T$$

$$\tau_b = 10^4 \log(e) \tau$$


Temperature correction

$$R_6 = \sum_{i=1}^n w_i F_c(\lambda_i) = \sum_{i=1}^n w_i F(\lambda_i) + \sum_{i=1}^n w_i \tau_b(\lambda_i) T$$

$$\sum_{i=1}^n w_i = 0$$

$$\tau'_b(\lambda_i) = \tau_b(\lambda_i) - \tau_b(\lambda_0)$$

Statistics from The EUBREWNET DB




Temperature characterisation of Brewer determined in the laboratory



GOBIERNO DE ESPAÑA



MINISTERIO DE AGRICULTURA, ALIMENTACIÓN Y MEDIO AMBIENTE





Description of the experiment


Measurements at the PTB


TH – Temperature and humidity sensor
SIC – Silicon carbide photodiode
CVC – Current-to-voltage converter
CC – Climate chamber
DVM – digital voltmeter



Measurements at the PTB



















Mode	Input Port	Duration (min)	Cycles	Filter1	Filter2
ozone	Internal Lamp	10	20	256	0
aod	Internal Lamp	20	10	256	0/64
uv	Internal Lamp	15	4	256	0
ozone	<u>Global Port</u>	5	50	128	0
aod	<u>Global Port</u>	15	10	128	0
uv	<u>Global Port</u>	10	1	128	0
ozone	Direct Port	5	20	256	0/64
aod	Direct Port	15	10	256	0/64
uv	Direct Port	10	1	256	64

Temperature characterisation of Brewer determined in the laboratory



$$R_6 = \sum_{i=1}^n w_i F(\lambda_i) \quad (\text{no temperature correction applied})$$



Temperature correction

$$F(\lambda_i) = F_c(\lambda_i) - \tau_b(\lambda_i)T$$

$$F(\lambda_i) - F(\lambda_0) = F(\lambda_i) - F(\lambda_0) - (\tau_b(\lambda_i) - \tau_b(\lambda_0))T$$

$$F(\lambda_i) - F(\lambda_0) = F(\lambda_i) - F(\lambda_0) - \tau'_b(\lambda_i)T$$







$$\tau'_b(\lambda_i) = \tau_b(\lambda_i) - \tau_b(\lambda_0)$$

Temperature characterisation of Brewer determined in the laboratory









GOBIERNO DE ESPAÑA







MINISTERIO DE AGRICULTURA, ALIMENTACIÓN Y MEDIO AMBIENTE




Temperature characterisation of Brewer determined in the laboratory




Temperature characterisation of Brewer determined in the laboratory







Comparison of Temperature coefficients obtained from laboratory measurements and field measurements




RBCC-E	Temperatures: 19 - 35	Dates: 2015-05-23 / 2015-06-02
Jun-Nov	Temperatures: 17 - 34	Dates: 2015-06-12 / 2015-11-30
Jul	Temperatures: 18 - 34	Dates: 2015-07-01 / 2015-07-31

Temperature characterisation of Brewer determined in the laboratory



Effect of PTFE diffuser on the instrument temperature dependence.









Ylianttila and Schreder [2005] "Temperature effects of PTFE diffusers"

Temperature characterisation of Brewer determined in the laboratory






GOBIERNO DE ESPAÑA

MINISTERIO DE AGRICULTURA, ALIMENTACIÓN Y MEDIO AMBIENTE



Temperature characterisation of Brewer determined in the laboratory





Conclusions

- The temperature sensitivity of Brewer spectrophotometer determined from internal lamp measurements and from external lamp measurements through the global port are in close agreement. Slits 5 and 6 present a non-linear dependence for temperature above 30°C.
- However, the temperature sensitivity determined from external lamp measurements through the direct port presents a linear behavior. This discrepancy with the other measurements may indicate a temperature dependence of the quartz window.
- Brewer global measurements are affected by a transmittance increase at about 20°C due to the PTFE diffuser. This effect is wavelength dependent, changing from 5% at 300nm to 3% at 360nm.

Thank you for your attention